

Prompt:

Create clear and informative instructional materials for [INSERT SUBJECT]. The instructions should include the following:

A clear goal of the instruction, which is the desired outcome for the user.

Breakdown the instructions into small, manageable steps using headings, bullet points, and numbered lists.

Use simple and clear language that is easy to understand, avoiding jargon or technical terms.

Include examples and visual aids such as images, diagrams, and videos.

Test the instructions with a small group of users and take feedback to make revisions as needed.

Review and proofread the instructions multiple times to ensure there are no errors or inconsistencies.

Make sure the instructions are in a format that is easy to access, like a printed manual, a webpage, or an app.

Keep the instructions up-to-date with current information and technology.

Begin your instructions with a message thanking the customer for buying the product.

ChatGPT Response:



Dear Flowflex COVID-19 Antigen Home Test user,

Thank you for purchasing our product. Our goal is to provide you with a convenient and reliable way to test for COVID-19 in the comfort of your own home. In order to achieve this goal, please follow the instructions below carefully.

Step 1: Preparing for the test

- Wash your hands thoroughly with soap and water.
- Open the test kit and remove the contents.
- Find a clean and flat surface to perform the test.
- Make sure the sample collection area is clean and dry.

Step 2: Collecting the sample

- Remove the swab from its packaging and hold it by the handle.
- Tilt your head back slightly and insert the swab into your nostril.
- Gently rotate the swab around for 10-15 seconds.
- Remove the swab from your nostril and repeat the process with the same swab in your other nostril.

Step 3: Testing the sample

- Open the test cartridge and insert the swab into the sample well.
- Press down on the swab handle to release the sample.
- Wait for 15 minutes and do not move or touch the test cartridge.
- After 15 minutes, read the test result.

Step 4: Reading the result

- Look at the test result window.
- If there is one line in the "C" region, the test is negative.
- If there are two lines in the "C" and "T" regions, the test is positive.
- If there is no line in the "C" region, the test is invalid and needs to be repeated.

Step 5: Disposing of the test

- Wrap the used swab in tissue or paper and dispose of it in the trash.
- Dispose of the test cartridge and packaging in the trash.

We hope these instructions are helpful and easy to follow. If you have any questions or concerns, please refer to our website or contact customer service for assistance.

Please note that the accuracy of the test may be affected by a variety of factors such as the timing of the test, sample collection and handling, and other factors that may impact the performance of the test.

Thank you for using our Flowflex COVID-19 Antigen Home Test. Stay safe and healthy!